

Multi-purpose Phantom Cradle

- Quicker phantom setup and positioning
 - Ideal for PET calibration phantoms
- Supports the ACR CT accreditation phantom
- Tape no longer needed to secure phantom



Phantom Cradle (NMTK-04)

This is a multi-purpose cradle that supports cylindrical QC phantoms with diameters of 16-24 cm. The cradle is only used for table top scanning. Cannot be used in a head support device or an MR head coil. Phantom length is not restricted due to the cradle design. Cradle is constructed from closed cell polyethylene foam. Dimensions: 30.5 cm (l) x 29 cm (w) x 9 cm (h). For SPECT scanning when radius is a consideration cradle NMTK-03 is recommended in place of NMTK-04.

